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Reference

1. Rodrick J. Treatment of Post-Surgical Edema in the Orthopaedic Patient – A case report. www.lymphnet.org/vol 22/2. Accessed on 21st December 2017. **2.** Sugrue M et al. World Union of Wound Healing Societies (WUWHS) Consensus Document. Closed surgical incision management: understanding the role of NPWT. Wounds International, 2016 **3.** Villeco JP. Edema: a silent but important factor. Journal of Hand Therapy. 2012;25(2):153-162. **4.** Ebert JR, Joss B, Jardine B, Wood DJ. Randomized trial investigating the efficacy of manual lymphatic drainage to improve early outcome after total knee arthroplasty. Archives of physical medicine and rehabilitation. 2013;94(11):2103-2111. **5.** Patel VP, Walsh M, Sehgal B, Preston C, DeWal H, Di Cesare PE. Factors associated with prolonged wound drainage after primary total hip and knee arthroplasty. The Journal of bone and joint surgery American volume. 2007;89(1):33-38. **6.** Saleh K, Olson M, Resig S, et al. Predictors of wound infection in hip and knee joint replacement: results from a 20 year surveillance program. Journal of Orthopaedic Research. 2002;20(3):506-515. **7.** Pua Y-H. The Time Course of Knee Swelling Post Total Knee Arthroplasty and Its Associations with Quadriceps Strength and Gait Speed. The Journal of Arthroplasty. 2015;30(7):1215-1219. **8.** Holm B, Kristensen MT, Bencke J, Husted H, Kehlet H, Bandholm T. Loss of Knee-Extension Strength Is Related to Knee Swelling After Total Knee Arthroplasty. Archives of Physical Medicine and Rehabilitation. 2010;91(11):1770-1776. **9.** Deakin AH, Iyayi-Igbinovia A, Love GJ. A comparison of outcomes in morbidly obese, obese and non-obese patients undergoing primary total knee and total hip arthroplasty. the surgeon. 2017 **10.** Chesney D, Sales J, Elton R, Brenkel IJ. Infection after knee arthroplasty: a prospective study of 1509 cases. The Journal of arthroplasty. 2008;23(3):355-359. **11.** Badia JM, Casey AL, Petrosillo N, Hudson PM, Mitchell SA. Impact of surgical site infection on healthcare costs and patient outcomes: a systematic review in six countries. Journal of hospital infection. Pages 1 -15. 2017. **12.** Karlakki SL, Hamad, A. K., Whittall, C., Graham, N. M., Banerjee, R. D., & Kuiper, J. H. Incisional negative pressure wound therapy dressings (iNPWTd) in routine primary hip and knee arthroplasties. Bone and Joint Research. 2016;5(8):328-337 **13.** Webb LX. New Techniques in Wound Management: Vacuum-Assisted Wound Closure. Journal of the American Academy of Orthopaedic Surgeons. 2002;10(5):303-311 **14.** Galat DD, McGovern SC, Hanssen AD, Larson DR, Harrington JR, Clarke HD. Early return to surgery for evacuation of a postoperative hematoma after primary total knee arthroplasty. The Journal of bone and joint surgery American volume. 2008;90(11):2331-2336. **15.** Nherera LM et al. Cost-effectiveness analysis of single-use negative pressure wound therapy dressings (sNPWT) to reduce surgical site complications (SSC) in routine primary hip and knee replacements. Wound Repair & Regeneration 2017;25(3):474-482. DOI:10.1111/wrr.12530. **16.** Boon NA, Colledge NR, Walker BR, Hunter JAA. Davidson's Principles & Practice of Medicine, 20th edition, publisher: Elsevier, published 2006. **17.** Chan PK, Brenkel IJ, Aderinto J. The outcome of total hip and knee arthroplasty in diabetics. Current Orthopaedics. 2005;19(1):59-67. **18.** Lumb, H.Data on file report.1102010. April 2011 **19.** Loveluck et al. Biomechanical modelling of the forces applied to closed incisions during single-use NPWT. Eplasty 2016 **20.** Pellino G, Scuiadone G, Candilio G, Campitiello F, et al., Effects of a new pocket device for Negative Pressure Wound Therapy on surgical wounds of patients affected by Crohn's disease: a pilot trial. Surgical Innovation. 2013. 21; (2):204–212. **21.** Hyldig, N., Birke-Sorensen, H., Kruse, M., Vinter, C., Joergensen, J., Sorsensen, J., Mogensen, O., Lamont, R. and Bille, C. (2016). Meta-analysis of negative-pressure wound therapy for closed surgical incisions. British Journal of Surgery, 103(5), pp.477-486.**22.** Canonico S et al. Therapeutic possibilities of portable NPWT. Initial mulitidisciplinary observation with the negative pressure device. Acta Vulnov. Vol 10, Pages 57 -66.2012 **23.** Malmsjö M, Huddleston EM, Martin R. Biological Effects of a Disposable, Canisterless Negative Pressure Wound Therapy System; Eplasty 2014;14:e15Page 4: PICO Intro **24.** DS.16.352.R Measurement of NPWT device air leak tolerance – PICO vs PICO II. September 2016. **25.** DS.16.351.R Comparison of noise levels generated by the PICO and PICO II devices in operation. September 2016. **26.** WMP.11446.UEF/R3 Project Fairbanks Human Factors Summary Report Issue 5. G. Walker, May 2017. **27.** Hurd, T; Gregory, H; Jones, A; Brown, S; A multicenter in-market evaluation of Allevyn Gentle Border; Wounds UK 2009, Vol 5, No 3 **28.** Data on file reference DS/11/037/R2 – In-vitro wound model testing of PICO at a moderate exudate flow rate; Sarah Roberts, March 2011, (in-vitro) **29.** Data on file reference DS/11/057/R2 – In-vitro wound model testing of PICO at a low exudate flow rate; Sarah Roberts, April 2011, (in-vitro) **30.** Data on file.DS/17/253/R version 2. Project Opal PICO 7 system stability testing, initial time point. November 2017 **31.** Data on file. DS/17/669/R. Negative pressure transmission through dressing port under a compressive force.November 2017 **32.** Galiano et al, The effects of a Single Use Canister-free Negative Pressure Wound Therapy (NPWT) sytem on the prevention of postsurgical wound complications in patients undergoing bilateral breast reduction surgery, 30th Annual Scientific Meeting, 2014 **33.** Strugala V, Martin R. Meta-analysis of comparative trials evaluating a prophylactic Single-use Negative Pressure wound Therapy system for the prevention of surgical site complications. Surgical Infections, Vol 18. 2017. **34.** Karlakki S et al., (2013) Negative pressure wound therapy for management of surgical incision in orthopaedic surgery, BJR, Vol 2, No. 12, 2013. **35.** Jenks PJ, et al. Clinical and economic burden of surgical site infection (SSI) and predicted financial consequences of elimination of SSI from an English hospital. J Hosp Infect 2014; 86: 24-33 **36.** Tanner J et al., Post discharge surveillance to identify colorectal surgical site infection rates and costs. Journal of hospital infection 72(3), PP. 242-250 (2009) **37.** Tanner 2011, Tanner et al., (2001), Surveillance of surgical sites in primary care, Nursing Times, Vol 107 No 9. **38.** Data on file report DS/14/065/R. Comparison of pressure transmission through port area of current PICO and PICO Soft Port designs. April 2014. **39.** Data on file report DS/14/066/R. Comparison of pressure transmission through RENASYS drape and PICO Soft Port dressings. May 2014. **40.** as Gristina AG.; Biomaterial-centered infection: microbial adhesion versus tissue integration . Science. 1987;237(4822):1588–95.

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©Duben 2018 Smith & Nephew 12356
Medical Ltd